

Enteral Feeding Bags

Available in both Gravity and Pump Bags

Pump Enteral Feeding Bags

Model iEPB 1000, 1000mL, 30 pieces per case

Gravity Enteral Feeding Bags

Model iEGB 1000, 1000mL, 30 pieces per case

- 1000-mL feeding bags
- Designed to be attached to a feeding pump or for gravity feeding
- Large, easy-open cap with locking tab for filling of Enteral formula
- Integral hanger for suspension of the bag
- Graduated capacity markings for accurate filling and measuring
- Extra long clear tubing with roller clamp for patient control of feeding rate
- Drip chamber for easy viewing of the liquid nutrition
- Non IV compatible Slip-Fit and ENFit connectors, with caution label attached
- Instructions for use are printed on the bag
- Pump bags are compatible with many popular pumps for Enteral feeding
- Non-sterile; ready to fill, hang and use
- Latex Free

The model EPB 1000 Pump Feeding Bag fits the following feeding brands of pumps:

- 1) ELAN PHARMA® Enteral Pumps: EP85, EP80, EP70, EP60
- 2) COMPAT® Enteral Pumps: 235, 225, 205
- 3) KANGAROO® Enteral Pumps: 330, 220, 222
- 4) OBRIEN KMI® Enteral Pumps: KM 70, KM 60, KM 50
- 5) Sentinel
- 6) Nestle/Zevex 2200
- 7) Medline® EntraFlo
- 8) Other popular pumps

