

IC-1 Disinfectant Cap



The IC-1 Disinfection Cap **DISINFECTS** all Female Needle Free Valves in less than one minute and provides bacterial **PROTECTION** for up to 7 days.



| | IC-1 | Competitor B | Competitor C |
|---------------------------|------|--------------|--------------|
| Disinfects in 1 minute | ✓ | ✓ | ✓ |
| Protects for up to 7 days | ✓ | ✓ | ✓ |
| 70% Isopropyl Alcohol | ✓ | ✓ | ✓ |
| 360* Rotation | ✓ | ✗ | ✓ |
| IV pole strips offered | ✓ | ✓ | ✗ |
| 10 caps per strip | ✓ | ✓ | ✗ |
| Individual caps offered | ✓ | ✓ | ✓ |

PACKAGED IN 2 CONFIGURATIONS


Model i-FIND-2000

Individually packed
 200 pieces to the box
 10 boxes to the case

Model i-FSTRIP-2500

On strips
 10 pieces on a strip
 25 strips in a box
 10 boxes to the case

Ergonomic Design for Universal Use



Proprietary 70% IPA
 Applicator