


SUCTION CANISTER

FOR USE WITH SUCTION TUBING SET **REF** iST

REF iSC600

QTY - 1EA



The iMed 600mL or 1000mL Suction Canister is a non-powered, single patient use device that consists of a manually operated plastic disposable evacuation system intended to provide a vacuum for suction drainage of body fluids.

- The canister is Manufactured from Shatterproof, Medical Grade Plastic for Safety.
- The Drainage Bottle has a Universal, Needle-Free Luer Lock Connection.
- Bottle Features Tapered Sides for Easy Handling
- Large Green Vacuum Indicator.
- Adjustable, Removable Hanging Strap



Read instructions prior to use



Single use



Do not use if pack is opened or damaged



Store in cool, dry area

Sterility guaranteed unless package is opened or damaged.

Not made with Natural Rubber Latex

LOT

